




SEDSIM[®]

The first true EEG Simulator for Anesthesia.



Module

- Portable
- 4 channel EEG
- Connection to MOC4



Head

- Compatible with different brands
- Connection to sensor



HOW TO USE IT




ARCOR
MEDICAL




-OR-



-OR-




Sedsim.exe



STATIC SIGNALS



ARCOR
MEDICAL



- EEG to generate stable BIS index
- Switch between different status
- Generate EEG bilateral or unilateral suppressions
- Simulate suppressions by low NIRS
- Change amplitude to simulate aged brains




USE CASES

- **Latency Stress Test:** Rapidly switch EEG states (e.g., suppression to awake) to measure response delays in BIS, SEF, SR, and suppression time across different monitors.
- **Hidden Screen Workshop:** Let users make therapeutic decisions (e.g., induction bolus) without seeing the monitor, then reveal EEG and index trends to assess prediction accuracy.
- **Spectrogram Amplitude Drill:** Compare EEG signals of equal frequency but different amplitudes to show how aging or low-voltage brains affect spectral displays and BIS readings.
- **Unilateral vs Bilateral Suppression:** Simulate asymmetric suppression to highlight differences in monitor interpretation and relevance for localized cerebral injuries.
- **Cerebral Hypoxia Simulation:** Introduce suppression due to low cerebral oxygenation (not drug effect) to demonstrate how indices may not reflect true clinical cause.
- **Index Drift Demo:** Provide a constant EEG signal and observe how BIS, SR, and SEF fluctuate over time, exposing monitor noise, filtering, and stability limitations.

CUSTOM SIGNALS



ARCOR
MEDICAL



- Create your own Delta, Theta, Alpha and Beta waves
- Customize your spectrogram
- Recreate different drug signatures
- Recreate spindles for intermittent states
- Recreate seizures


USE CASES

- **Drug Specific EEG Signatures:** Simulate characteristic EEG patterns produced by common agents such as propofol, volatile anesthetics, ketamine, and dexmedetomidine to train recognition. **Index Sensitivity to Waveform Changes:** Explore how changes in specific EEG waveforms, including alpha, delta, and theta, affect BIS and related depth of anesthesia indices. **Alpha Dropout Variants:** Demonstrate different types of alpha wave loss such as progressive thinning, amplitude reduction, or abrupt disappearance. **Ketamine and Propofol Interaction:** Reproduce EEG changes seen after a ketamine bolus during continuous propofol infusion, including spectrum widening, enhanced beta activity, and altered alpha band patterns. **Dexmedetomidine EEG Patterns:** Simulate low dose dexmedetomidine effects featuring intermittent spindle-like alpha bursts followed by dominant delta activity to illustrate microstate transitions during sedation. **Seizure Pattern Simulation:** Introduce bursts of irregular mixed frequency activity mimicking seizure onset, followed by a return to normal EEG, to practice identifying abnormal events within anesthetic EEG recordings.
-

CLINICAL CASES



ARCOR
MEDICAL



- Ideal for booths or unattended monitors on display
- 20 minute cases on repeat
- 3 representative cases:
 - Typical TIVA anesthesia case
 - Ketamine to Propofol transition
 - Unilateral suppressions on left hemisphere from
- Load EDF case: Recreate a pre-recorded real case from the OR


CLINICAL CASE 2



ARCOR
MEDICAL



- Patient with Ketamine infusion
 - Delta and slight Beta waves
- Propofol bolus
 - Alphawaves appear
- 2nd Propofol bolus
 - Suppressions appear



ADVANCED MODE



ARCOR
MEDICAL



STATIC SIGNALS | CLINICAL CASES | CUSTOM SIGNALS | **ADVANCE MODE**

PATIENT INFORMATION

Gender: MALE

Weight (Kg): 70

Height (cm): 170

Age: 30

ANESTHETIC INFORMATION

Drug: PROPOFOL

Infusion mode: TCI

Model: SCHNIDER

Ce/Cp/Ca: EFFECT_SITE

Target ($\mu\text{g}/\text{mL}$): 0.0

Induction time (min): 1.0

Maintenance time (min): 1.0

Dilution (mg/ml): 10.0

Initial speed (mL/h): 0.0

Initial bolus (ml): 0.0

START **STOP**

Time: 00:00:30

Information

SIGNAL LOADED...

Suppression: 0%

Configuration

Amplitude: 100%

NIRS: 80%

DSA Lower: -40

DSA Upper: 40

Delta: 38.7% Theta: 10.0% Alpha: 27.8%
Beta: 19.9% Gamma: 3.5%

- Simulate pharmacokinetics and their effects on EEG
- Effect site and plasma concentration simulation
- PK/PD models for propofol, volatile, dexmedetomidine and Ketamine
- Bolus, continuous infusion and TCI modelling for infusion modes

STATIC SIGNALS | CLINICAL CASES | CUSTOM SIGNALS | **ADVANCE MODE**

PATIENT INFORMATION

Age: 30 years
Gender: MALE
Height: 170 cm
Weight: 70 Kg

ANESTHETIC INFORMATION

Anesthetic: PROPOFOL
Dilution: 10.0mg/ml
Infusion mode: TCI
Model: SCHNIDER
Ce/Cp: EFFECT_SITE

SIMULATION INFORMATION

Target: 4.0
Ce[]: 4.0
Cp[]: 4.1
Inf Rate (ml/h): 73

SEND/CHANGE

New Target ($\mu\text{g}/\text{mL}$): 4.0

New speed (mL/h): 0.0

bolus (ml): 0.0

START **STOP**

Time: 00:05:59

Information

ACTIVE SIMULATION...

Suppression: 0%

Configuration

Amplitude: 100%

NIRS: 80%

DSA Lower: -40

DSA Upper: 40

Delta: 84.0% Theta: 5.5% Alpha: 10.1%
Beta: 0.4% Gamma: 0.0%



USE CASES

Compare Infusion Strategies in Real Time: Simulate the EEG response differences between bolus, manual infusion, and TCI-based administration for the same drug, helping clinicians visualize how delivery mode affects anesthetic depth and EEG dynamics. **Effect Site vs Plasma Lag Demonstration:** Show the temporal disconnect between plasma and effect-site concentrations for agents like propofol, using EEG to illustrate why immediate EEG changes do not always match dosing changes. **Personalized Dosing Scenarios:** Input patient-specific parameters (e.g., age, weight, sex) into PK models and observe differences in drug onset, peak, and EEG response. Use this to teach inter-patient variability and the rationale for closed-loop systems. **Overdose and Underdose Training:** Deliberately simulate excessive or insufficient dosing scenarios (e.g., overly rapid propofol TCI ramp) to show EEG suppression or inadequate anesthesia, training users to recognize red flags and correct dosing plans. **Emergence Timing and Planning:** Simulate different weaning strategies from TCI or volatile agents and compare EEG trajectories toward wakefulness, emphasizing how PK/PD models can guide smoother emergence without awareness or agitation.



SEDSIM

4-channel EEG
Simulator for
Anesthesia.



SEDSIM

4-channel EEG
Simulator for
Anesthesia.



ARCOR

MEDICAL

